

Vacuum control unit VC 900



Technical features

- Control of the vacuum application
- Separate control unit with pressure sensors and two-step controlled valve to be placed independently from the operating unit
- Easy to use

Technical data

Operating unit

Housing materials	Plastic
Power supply	Power cord, cable length 150 cm
Protection class	IP30
Mains voltage [V]	100 – 240
Frequency [Hz]	50/60
Operating current [A]	max. 1.0
Dimensions W x D x H [mm]	181 x 101 x 67
Accuracy of measurement	± 1 mbar
Measuring limit	upper: 1100 mbar abs. lower: 0 mbar abs.
Interface	Mini-USB

Control unit

Dimensions W x D x H [mm]	155 x 60 x 109
Hose connections:	
- for gas in	ID10, PVDF
- for gas out	ID10, PVDF
- for venting	ID4, nickel-plated brass

Other data

Total weight [lbs]	2.6
Permissible media and ambient temperature [°C]	+10 to +40

Installation

Operating unit can be used as bench-top device providing 3 different display positions

Operating unit can be mounted at the wall

Dimensions [mm]

